

INSTRUCTIONS FOR USE



■ BD MAX™ CT/GC/TV 20-Day QC Panel

INTENDED USE

The BD MAX™ CT/GC/TV 20-Day QC Panel is intended for use as an external assayed positive quality control material to monitor the performance of *in vitro* laboratory nucleic acid testing procedures for the qualitative detection of *Chlamydia trachomatis*, *Neisseria gonorrhoeae* and *Trichomonas vaginalis* with the BD MAX™ CT/GC/TV Assay on the BD MAX™ System. The controls comprise cultured and inactivated *C. trachomatis*, *N. gonorrhoeae* and *T. vaginalis*.

The BD MAX CT/GC/TV™ 20-Day QC Panel is not intended to replace manufacturer controls provided with the device.

SUMMARY AND EXPLANATION

The BD MAX™ CT/GC/TV 20-Day QC Panel is an easy-to-use external quality control that can be used to monitor the extraction, amplification and detection of the BD MAX™ CT/GC/TV assay. Routine use of quality controls helps in identifying problems with the test system due to instrument or reagent failure. At a minimum, controls should be run at the frequency recommended by the assay system manufacturer.

COMPOSITION

Each BD MAX™ CT/GC/TV 20-Day QC Panel consists of 20 individually packaged positive control pellets. Each individually packaged pellet consists of inactivated *Chlamydia trachomatis*, *Neisseria gonorrhoeae* and *Trichomonas vaginalis* pathogens at target levels designed to provide reproducible performance above the limit of detection of the assay. The organisms are prepared in a buffered suspension with materials of animal origin, preservatives and stabilizers. The suspension is lyophilized into a ready-to-use pellet.

WARNINGS AND PRECAUTIONS

- For In Vitro Diagnostic use only.
- Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner.
- Do not open foil pouch until ready to use.
- This product must be treated as a viable specimen and handled in accordance with Biosafety Level 2 practices as described in the United States Department of Health and Human Services Centers for Disease Control and Prevention (CDC) and National Institutes of Health (NIH), Biosafety in Microbiological and Biomedical Laboratories, or other equivalent guidelines.
- Wear proper personal protective equipment.
- Refer to the Safety Data Sheet (SDS) for more detailed information. The SDS can be located on the Microbiologics website at www.microbiologics.com or by contacting Technical Support at 320.229.7045 or U.S. Toll Free 1.866.286.6691.
- These products are not made with natural rubber latex.



STORAGE AND EXPIRATION

Store the Inactivated Helix Elite Molecular Standards at 2°C - 25°C in the original packaging up to the indicated expiration date. After opening the foil pouch, use the pellet within 5 hours.

Inactivated Helix Elite Molecular Standards should not be used if:

- Stored improperly
- There is evidence of excessive exposure to heat or moisture
- The packaging is damaged or shows evidence of tampering
- The expiration date has passed

INSTRUCTIONS FOR USE

Note: Follow the Specimen Preparation section of the instructions for use in the BD MAX™ CT/GC/TV assay package insert. Do not follow the Quality Control instructions for sample preparation. The Microbiologics controls are ready to use controls and must be processed as a patient sample after transfer to the UVE Sample Buffer Tube.

1. Read the instructions for use specific to the BD MAX™ CT/GC/TV assay.
2. Tear open the Helix Elite™ pouch at notch. Remove the tube from the pouch and ensure that the pellet is at the bottom of the vial before opening the tube.
3. Uncap the tube and carefully transfer the pellet to the UVE Sample Buffer Tube. **DO NOT DILUTE.**
4. Vortex for 5 seconds on high to mix.
5. Invert the UVE Sample Buffer Tube 5 times.
6. Vortex again for 5 seconds on high to mix.
7. Pre-warm and process following instructions listed in the BD MAX™ CT/GC/TV instructions for use.
Note: The UVE Sample Buffer Tube should only be tested once.
8. Mixture must be used within 5 hours of mixing.

MATERIALS REQUIRED BUT NOT PROVIDED

- BD MAX™ System
- BD Pre-warm Heater Kit
- BD MAX™ UVE Sample Buffer Tube from the UVE Specimen Collection Kit
- BD MAX™ CT/GC/TV assay
- BD MAX™ PCR Cartridges
- Package insert of the BD MAX™ assay
- Vortex mixer

LIMITATIONS

- Quality control materials should be used in accordance with local, state, federal regulations, and accreditation requirements.
- Once inoculated with the control, the UVE Sample Buffer Tube should only be processed once on the BD MAX™ System. Studies have not been performed to support performance of the control if the remaining suspension in the UVE Sample Buffer Tube is stored and processed a second time.

EXPECTED VALUES

The following results are expected when running the BD MAX™ CT/GC/TV 20-Day QC Panel with the BD MAX™ CT/GC/TV assay:

Expected Assay Result	Interpretation
CT POS	<i>Chlamydia trachomatis</i> DNA detected.
GC POS	<i>Neisseria gonorrhoeae</i> DNA detected.
TV POS	<i>Trichomonas vaginalis</i> DNA detected.

A single test of the BD MAX™ CT/GC/TV 20-Day QC Panel should yield the expected results when tested. A positive control that yields a negative result may be indicative of a specimen handling/preparation problem. Repeat the test using a new pellet and a new UVE Sample Buffer Tube.

The performance of the BD MAX™ CT/GC/TV 20-Day QC Panel was evaluated in a study that was performed at three testing sites using three different production lots on three different BD MAX™ Systems. A total of six operators performed the testing. The results of the study are summarized below.

Analyte	Agreement (%) by Test Site/BD MAX System			
	Site 1 ¹	Site 2 ¹	Site 3	Overall
<i>C. trachomatis</i>	30/30 (100)	30/30 (100)	30/30 (100)	90/90 (100)
<i>N. gonorrhoeae</i>	30/30 (100)	30/30 (100)	29/30 (97)	89/90 (99)
<i>T. vaginalis</i>	30/30 (100)	30/30 (100)	30/30 (100)	90/90 (100)

¹Three Unresolved results were obtained; in all cases a new control was retested and the expected results were obtained

KEY OF SYMBOLS



Batch Code (Lot)



Catalog Number



Caution, Consult Accompanying Documents



Contains sufficient for <n> tests



In Vitro Medical Device



Manufacturer



Temperature Limitation



Use By



Refer to Instructions for Use



Telephone Number

PRODUCT WARRANTY

- These products are warranted to meet the specifications and performance printed and illustrated in product inserts, instructions, and supportive literature.
- The warranty, expressed or implied, is limited when:
 - The procedures employed in the laboratory are contrary to printed and illustrated directions and instructions

- The products are employed for applications other than the intended use cited in product inserts, instructions, and supportive literature

NOTICE TO PURCHASERS _____

The purchase of this product allows the purchaser to use it for Research and Quality Control. No general patents or other license of any kind other than this specific right of use from purchase is granted hereby. No other rights are conveyed expressly, by implication or by estoppel to any other patents. Furthermore, no rights for resale are conferred with the purchase of this product.

The Microbiologics logo and Helix Elite are registered trademarks of Microbiologics, Inc.

WEBSITE _____

Visit our website, www.microbiologics.com, for current technical information and product availability.

BIBLIOGRAPHY _____

1. BD MAX™ CT/GC/TV [package insert]. Becton, Dickinson and Company, Sparks, MD 21152 USA.

ACKNOWLEDGEMENTS _____



Microbiologics, Inc.
200 Cooper Avenue North
St. Cloud, MN 56303 USA
www.microbiologics.com

Customer Service
Tel. 320-253-1640
U.S. Toll Free 800-599-BUGS (2847)
Email: info@microbiologics.com

Technical Support
Tel. 320-229-7045
U.S. Toll Free 866-286-6691
Email: techsupport@microbiologics.com

ILLUSTRATED INSTRUCTIONS

1

Read the instructions for use specific to the BD MAX™ CT/GC/TV assay.


2

Tear open the Helix Elite™ pouch at notch. Remove the tube from the pouch and ensure the pellet is at the bottom of the vial before opening the tube.



3

Uncap the tube and carefully transfer the pellet to the UVE Sample Buffer Tube.




4

Vortex for 5 seconds on high to mix.




5



Invert the UVE Sample Buffer Tube 5 times.

6

Vortex again for 5 seconds on high to mix.



7

Pre-warm and process following instructions listed in the BD MAX™ CT/GC/TV instructions for use.

Note: The UVE Sample Buffer Tube should only be tested once.

8

Mixture must be used within 5 hours of mixing.