

## Preservative Testing – Choice of Challenge Isolates

By Dr. Philip Geis



**"What you see is that the most outstanding feature of life's history is a constant domination by bacteria."~Stephen Jay Gould**

An effective preservative capability is a central element of drug, cosmetic and household and institutional product microbiological risk assessment and quality. Such a capability can mitigate both incidental low-level contamination from a controlled manufacturing system, as well as introduction of contamination through typical, expected product use; establishing appropriate product quality for the life of the product. Whereas some products do not require frank preservation by design (e.g. aseptic production/single use), formulation (e.g. high ethanol content, extremes of pH and anhydrous) or constraining dynamics of use (e.g. near term expiration dates, refrigeration), chemical preservation and preservative qualification remain a primary element of quality for the majority of aqueous consumer products.

### **Preservative Effectiveness Testing**

As is seen with cosmetics, a relatively small set of preservatives are used to establish and maintain microbiological quality (1). Unfortunately simple addition alone or in combination is not sufficient to assure efficacy. The critical preservation capability must be demonstrated by preservative effectiveness test (PET) protocols that, as with preservatives themselves, is product category specific. Relevant PET's are USP <51> for drugs (2), AOAC protocol developed by the Personal Care Products Council for cosmetics (3), and a variety of methods (e.g. 4,5) for household and institutional products. A central element of each tests is consistency driving to a level of preservative efficacy that accommodates realities of microbial contamination risk in manufacturing

and product use.

These tests are relatively similar in that each exposes the product to an exaggerated microbial challenge, monitors survival over a 4 week period and derives a pass/fail determination-based comparison of microbial reduction observations to established criteria. The primary functional difference between category methods is the inoculum - the microbes used.

### Drug PET

As drug products are typically made under conditions of stringent hygiene and enjoy protective packaging, storage and expiration labeling, USP <51> with its small number of clinical isolates and modest kill rate expectations is more aligned to low-level potential contamination in making and use. The isolates indicated by <51> (Table 1) are primarily of clinical origin and sourced from the American Type Culture Collection. With original isolations ranging from 50 to almost 80 years ago, these isolates can be considered laboratory-adapted strains. None apparently bears resistance plasmids (6) but little additional information is available. The combination of GMP's, USP 51-driven efficacy, packaging and typical use appears to have been successful in maintaining a high standard of overall microbiological quality in for this product category.

**Table 1: USP PET Challenge Microorganisms**

Isolate	ATCC® No.	Date Submitted	Source	Comments
<i>Staphylococcus aureus</i>	6538™*	1930/1940	Clinical FDA	<ul style="list-style-type: none"> <li>• BSL 2</li> <li>• PEG's are inhibitory</li> <li>• Some potential for resistance to high pH</li> <li>• Min. Aw for growth 0.86</li> <li>• Agar grown cells more sensitive than broth grown (cationic biocides)</li> </ul>
<i>Pseudomonas aeruginosa</i>	9027™*	1943	Clinical	<ul style="list-style-type: none"> <li>• BSL 2</li> <li>• Similar disinfectant resistance to hospital isolates</li> <li>• Agar grown cells more sensitive than broth grown (cationic)</li> <li>• Broth grown cells give more reproducibility in preservative testing</li> </ul>

<i>Escherichia coli</i>	8739™*	1949	Clinical	<ul style="list-style-type: none"> <li>• BSL 1</li> <li>• Centrifuged/washed cells more sensitive to PHMB</li> <li>• Agar grown cells more sensitive than broth grown (cationic)</li> <li>• Broth grown cells more reproducible in preservative tests</li> </ul>
<i>Candida albicans</i>	10231™*	1950	Clinical	<ul style="list-style-type: none"> <li>• BSL 2</li> <li>• &gt; UV resistance than USP bacteria</li> <li>• Mycelial formation at acid pH</li> </ul>
<i>Aspergillus brasiliensis</i>	16404™*	1965	Botanical (blueberry)	<ul style="list-style-type: none"> <li>• BSL 1</li> <li>• Potential for biocorrosion</li> </ul>

### Cosmetic PET

Although some have used the USP method for qualification of cosmetic preservative systems (7), most cosmetic manufacturers use protocols based on the compendial method developed and qualified by the Microbiology Committee of the Personal Care Products Council (3). This protocol uses the USP 51 isolates and additional microbial isolates representative of species reported as manufacturing or in-use contaminants. A recent casual poll of major manufacturers found most using modified PCPC methodology with in-house manufacturing and consumer return isolates such as those described by Brannan et al. (8). This is not a trivial consideration as clinical isolates may be quite different from those recovered from the general environment (9). Unfortunately, commercial availability of authentic isolates is currently very limited, so most of these originated from clinical origins as with the USP isolates (Table 2). A rationale behind such additions would be that these microbes may persist at low levels in manufacturing systems so specific practical efficacy is assured. Further, some wish to establish of efficacy in context of isolates capable of developing tolerance or even resistance to preservatives. Though such capabilities may not be immediately expressed in subsequent laboratory culture, the genetic capability is presumably sustained.

Overall, cosmetic microbiology appears to have established the greatest control over microbial contamination. Reviews of FDA recall records for microbial contamination (10,11) have consistently found that annual recalls involving cosmetics were substantially fewer than those for foods or drugs.

Table 2: PCPC Additional PET Challenge Microorganisms

Isolate	ATCC® No.	Date Submitted	Source	Comments
<i>Enterobacter gergoviae</i>	33028™*	1948	Clinical (urine) CDC	<ul style="list-style-type: none"> <li>• BSL 2</li> <li>• Type species</li> <li>• France - Institute Pasteur</li> </ul>
<i>Burkholderia cepacia</i>	25416™*	1969	Botanical (onion)	<ul style="list-style-type: none"> <li>• BSL 2</li> <li>• Similar disinfectant resistance to hospital isolates</li> <li>• Cepobactin production (nonclinical isolate)</li> <li>• Lipase production</li> </ul>
<i>Acinetobacter baumannii</i>	19606™*	1966	Clinical (urine)	<ul style="list-style-type: none"> <li>• BSL2</li> <li>• Type species</li> <li>• Poor survival on dry environmental surfaces</li> <li>• Biofilm forming</li> </ul>
<i>Klebsiella pneumoniae</i>	10031™*	1946	Clinical FDA	<ul style="list-style-type: none"> <li>• BSL 2</li> <li>• Aw 0.96 is limit for growth,</li> <li>• Plasmolysis observed at Aw 0.93, rapid loss of viability at 0.88</li> <li>• Disinfectant resistance similar to hospital isolates</li> <li>• Antibiotics assay EP/USP/AOAC</li> <li>• Poor to no capsule formation</li> </ul>
<i>Staphylococcus epidermidis</i>	12228™*	1955	Clinical FDA	<ul style="list-style-type: none"> <li>• BSL 1</li> <li>• Plasmids</li> <li>• Lipase and chitinase production</li> <li>• Does not form biofilm</li> <li>• Heavy metal resistant</li> <li>• Antibiotics assay</li> </ul>

USP/EP/AOAC

### **PET for Everything Else**

As household and institutional products are rarely produced on manufacturing systems that establish effective microbiological control, the primary microbial component of challenge testing are manufacturing isolates. Base, unpreserved formulations of most household products would “pass” a USP 51 protocol due to formula characters including high or low pH, surfactant and salt concentrations, etc. So clinical isolates are less relevant to qualification of preservative efficacy for these products as their basic formulation, packaging and intended use militate against contamination in use. Contrast this with manufacturing where continuous challenge in absence of GMP’s that allows adaptation under relatively harsh conditions of microbes some might see as extremophiles (12,13,14). In-use contamination is a concern primary for those products whose use demands aqueous dilution. Even more than cosmetics, manufacturers in these categories rely on authentic manufacturing isolates to validate preservative efficacy.

### **Conclusion**

Successful risk assessment of any product requires consideration of manufacturing risk, preservation, package design and intended product use. Design of preservative efficacy testing, and most importantly selection of isolates against which efficacy is demonstrated, is central to establishing and maintaining product microbial quality. Compendial inocula drive to standard measure of general efficacy and many supplement such challenge with isolates more aligned to specific microbial risks relevant to formulation, manufacturing and product use. Though some attempts have been made to source such isolates (15), these remain primarily remain in-house and proprietary. Systems developed in absence of a clear understanding of microbial risk in making, use and anticipated misuse may well suffer substantial and surprising contamination.

## Microbiologics Custom Solutions



With Microbiologics Custom Solutions, our team of experts will professionally preserve your environmental isolate strain and create a product that will save your laboratory a tremendous amount of time, labor and money. Choose from our extensive catalog of product configurations, CFU concentrations and package sizes to design an easy-to-use, convenient solution that meets your exact needs.

Visit [www.microbiologics.com](http://www.microbiologics.com) to learn more about Microbiologics Custom Solutions.

## EZ-PEC™

Ready to use EZ-PEC™ quantitative microorganism preparations are designed to provide the required microorganism challenge concentrations for Antimicrobial Effectiveness Tests and Preservative Efficacy Tests. When processed as directed, EZ-PEC™ delivers a final concentration of  $1.0 \times 10^5$  to  $1.0 \times 10^6$  CFU per ml of the product being tested. For added convenience, a peel-off Certificate of Assay is included for quality control documentation.



## Biography



Dr. Geis earned a PhD in microbiology and mycology from the University of Texas. His career in microbiology began at a clinical lab in the US Army, moving to commercial media production, and in 1981 to The Procter & Gamble Company (P&G). Through almost three decades with P&G microbiology, Phil managed preservative and disinfectant development and studies of household and skin microbial ecologies and hygienic manufacturing. He was the first recipient of P&G's namesake award – Dr. Philip Geis Microbiology Quality Award. Dr. Geis brings unique global expertise and experience in diverse regulatory, manufacturing, product quality and consumer realities for a broad range of products from OTC drugs to fabric softeners to dog food.

## References

1. Geis PA. 2006. *Cosmetic Microbiology. A Practical Handbook.*
2. Sutton S, Porter D 2002. Development of the Antimicrobial Effectiveness Test as USP Chapter <51>. PDA J Pharma Sci Technol. 56:301-11.
3. Machtiger NA, et al. 2001. Determination of the Efficacy of Preservation of Non-Eye Area Water-Miscible Cosmetic and Toiletry Formulations: Collaborative Study. J AOAC Int. 84:101.
4. Geis PA, Rook A. 2011. Microbiological quality of consumer and industrial household products. Happi May/June: 82-7.
5. Cooke PK, et al, 1991. Preservative evaluation: designing an improved system. J Coat Technol 63: 33-8.
6. deSolis NMG, et al. 1994. Effect of plasmids conferring preservative resistance on performance of bacterial strains in compendial preservative efficacy tests. Euro J Pharma Sc. 2: 221-8
7. CTFA. 1990. CTFA survey: Test methods companies use. Cosm Toil 105:79-82.
8. Brannan DK, et al. 1990. Type of closure prevents microbial contamination of cosmetics during consumer use. Appl. Environ. Microbiol 56:1476-9.
9. Sutton S, Jimenez L. 2012. A review of reported recalls involving microbiological control 2004-2011 with emphasis on FDA considerations of "objectionable organisms." Am Pharma Rev 15:42-57.\
10. <http://www.fda.gov/Safety/Recalls/EnforcementReports/default.htm>
11. Beadle IR, et al. 1995. Studies on the Growth of Klebsiella sp. in a Cleaning Fluid. Int Biodeter Biodegrad 36:468-9.
12. Frenandez, P. Solvent tolerance in bacteria: role of efflux pumps and cross-resistance with antibiotics. Int J Antimic Agents 22:211-6
13. Nasser, R.M., et al. Outbreak of Burkholderia cepacia Bacteremia Traced to Contaminated Hospital Water Used for Dilution of an Alcohol Skin Antiseptic. Infect Cont Hosp Epidem 25: 231-9.
14. Decicco BT, et al. 1982. Factors affecting survival of Pseudomonas cepacia in decongestant nasal sprays containing thimerosal as preservative. J Pharm Sci. 71:1231-4.
15. Saxena, A.K., et al. 2013. MYCOsoft: A mycological database. J Bioinform Seq Anal 5:1-9.