


Microbiologics Dilutions Guide

This document outlines how to perform dilutions when using Microbiologics products. **Dilution is the process of making a solution weaker or less concentrated.** In microbiology, serial dilutions (log dilutions) are used to decrease a bacterial concentration to a required concentration for a specific test method, or to a concentration which is easier to count when plated to an agar plate. This document was created to provide a better understanding of dilutions and should be used as a guideline, not a replacement for laboratory procedures.

Log Dilutions

A log dilution is a tenfold dilution, meaning the concentration is decreased by a multiple of ten. To complete a tenfold dilution, the ratio must be 1:10. The 1 represents the amount of sample added. The 10 represents the total size of the final sample. For example, a sample size of 1 ml is added to 9 ml of diluent to equal a total of 10 ml.

Example: 1:10 dilution - if the concentration is 1,000 CFU, a one log dilution will drop the concentration to 100 CFU.



Decimal Numbers vs Scientific Notation

Decimal numbers can be converted to scientific notations by moving the decimal place the same number of places as the exponential number.


$$100.0 = 1.0 \times 10^2 \quad 1. \underset{\curvearrowright}{0} \underset{\curvearrowright}{0} = 100$$

$$350.0 = 3.5 \times 10^2 \quad 3. \underset{\curvearrowright}{5} \underset{\curvearrowright}{0} = 350$$

Decimal Number	1	10	100	1,000	10,000	100,000	1,000,000
Scientific Notation	10^0	10^1	10^2	10^3	10^4	10^5	10^6

Multiple Dilutions

Multiple dilutions are required to decrease the sample concentration by multiple logs. If the concentration is 35,000 CFU/ml (10^4), and 35 CFU/ml is the target concentration, the following serial dilutions can be performed.




Larger Dilutions

Decreasing the concentration using fewer dilutions is possible with the use of large volume dilutions. This can be done by performing a 1:100 dilution instead of 1:10. An example of this can be observed in the Epower™ instructions for membrane filtration using Microbiologics E3 Epower™ product.

The E3 Epower™ product provides 10^3 CFU per pellet which equates to 1 pellet in 1 ml equaling 10^3 CFU/ml. Placing a 10^3 Epower™ pellet in 10 ml will drop the concentration to 10^2 – this is a one log dilution.

A 1:100 dilution can be created by placing 1 pellet in 99 ml as instructed in the membrane filtration instructions. This will drop the concentration two logs from 10^3 to 10^1 CFU/ml.



EZ-CFU™ Dilution Example

When EZ-CFU™ is used according to directions, the following dilutions are conducted to reach a desired concentration of 10-100 CFU/0.1 ml.

1. Two pellets are placed in 2 ml of hydrating fluid = 1,000 - 10,000 CFU/ml.
2. 1:10 dilution is performed by placing 1 ml of the re-hydrated pellet solution into 9 ml of buffer = 100 - 1,000 CFU/ml.
3. 0.1 ml of the organism suspension plated to an agar = 10 - 100 CFU per 0.1 ml.