


# Cepheid GeneXpert® Controls



Microbiologics, the leading provider of control strains, offers QC panels designed specifically for Cepheid GeneXpert® systems. Our QC Panels feature individually packaged, pre-mixed pools of inactivated control strains recommended by the instrument manufacturer. The organisms are dried onto a swab to mimic the processing of a patient sample.

## Accurate, Reliable Results

- Recommended for routine QC
- Protects your reputation by providing documentation and assurance that the laboratory can consistently produce accurate results
- Improves daily workflow
- Ensures testing procedures and materials are working properly

## Easy-to-Use

- No mixing required
- Individually packaged and color coordinated labeling to simplify use and avoid cross contamination
- One simple catalog number for fast, easy ordering
- Room temperature storage
- For in vitro diagnostic use (IVD)

## INACTIVATED HELIX ELITE CONTROLS

*Each panel includes six Helix Elite™ swabs of each positive and negative control, and features color coordinated packaging for easy identification.*

**Cepheid Xpert® C. difficile/Epi Control Panel - Cat #8200**  
For use with the Xpert® C. difficile/Epi Assay

**Cepheid Xpert® CT/NG Control Panel - Cat #8188**  
For use with the Xpert® CT/NG Assay

**Cepheid Xpert® GBS LB Control Panel - Cat #8194**  
For use with the Xpert® GBS LB Assay

**Cepheid Xpert® MRSA NxG Control Panel - Cat #8195**  
For use with the Xpert® MRSA NxG Assay

**Cepheid Xpert® SA Nasal Complete Control Panel - Cat #8196**  
For use with the Xpert® SA Nasal Complete Assay

**Cepheid Xpert® Respiratory Control Panel - Cat #8199**

Contact us to learn more  
320.253.1640  
800.599.BUGS(2847)

[www.microbiologics.com](http://www.microbiologics.com)

 **Microbiologics®**  
A safer, healthier world.

FLY109 Rev H