

Burkholderia cepacia Now Available in EZ-Accu Shot™

Protect your customers from this notorious environmental contaminant



EZ-Accu Shot for fast, reliable Growth Promotion Testing

- Ready-to-use, no pre-incubation
- Instant dissolve pellet, no dilutions required
- Delivers 10-100 CFU per 0.1 ml as required by USP/ Ph. Eur./JP
- 50 tests per kit - ideal quantity for all labs
- Easy and economical refrigerated storage
- Up to 8 hours stability for ultimate flexibility
- Strains are four passages or fewer from reference culture
- Traceability to reference culture ensures authenticity
- Online Certificate of Analysis
- Product warranty

Available in 29 strains including:

- 0488A *Burkholderia cepacia* derived from ATCC® 25416™*
- 0392A *Aspergillus brasiliensis* derived from ATCC® 16404™*
- 0486A *Bacillus subtilis* subsp. *spizizenii* derived from ATCC® 6633™*
- 0443A *Candida albicans* derived from ATCC® 10231™*
- 0484A *Pseudomonas aeruginosa* derived from ATCC® 9027™*
- 0485A *Staphylococcus aureus* subsp. *aureus* derived from ATCC® 6538™*
- 0483A *Escherichia coli* derived from ATCC® 8739™*
- 0419A *Propionibacterium acnes* derived from ATCC® 11827™*

Visit microbiologics.com for a complete list of strains



Not intended for clinical use



*Look for the ATCC Licensed Derivative® Emblem for products derived from ATCC® cultures. The ATCC Licensed Derivative Emblem, the ATCC Licenses Derivative Word mark, and the ATCC Catalog Marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.