

May 2017

Dear Microbiologics Customer,

The name of Microbiologics product 0629 (derived from ATCC® 49148™\*) has been changed from *Staphylococcus xylosus* to *Staphylococcus gallinarum*. ATCC reported to Microbiologics that the name of ATCC® 49148™\* was changed after quality control testing showed 99% identification to *S. gallinarum* by BioMerieux VITEK GP, which was confirmed by 16S ribosomal gene sequencing (99.89% similarity). *S. gallinarum* was not in the API database at the time, so API testing showed closest match to *S. xylosus* at 95.6%.

Sincerely,

Jane Johnson  
Technical Support Specialist  
Microbiologics  
320.253.7400  
[jjohnson@microbiologics.com](mailto:jjohnson@microbiologics.com)



\*Look for the ATCC Licensed Derivative® Emblem for products derived from ATCC® cultures. The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.