



Extracting DNA or RNA for Molecular Testing?

Now offering a new universal lysis buffer for nucleic acid extraction



Easy-to-use UNEX Buffer* is a guanidium thiocyanate lysis buffer specially designed for universal nucleic acid extraction (UNEX) of DNA and RNA from viruses, bacteria, and parasites. Researchers from the Centers for Disease Control and Prevention developed and optimized the UNEX Buffer formulation to provide simple, scalable sample processing for molecular analysis.

Develop your own protocol

Versatile

- Optimized to work with a variety of pathogens and sample types

Scalable

- Compatible with column or phase-based extraction

Ready-to-Use

- No formulation or mixing required

Complex sample extraction

- Rotavirus
- *Cryptosporidium*
- *Giardia*
- Norovirus
- *Salmonella* spp
- Adenovirus
- And more!

* This product is not intended to diagnose, treat, cure or prevent disease.

Catalog # MR0501 UNEX Buffer

Bottle contains 100 ml of lysis buffer

www.microbiologics.com

Contact us to learn more
320.253.1640
800.599.BUGS(2847)

 **Microbiologics®**
A safer, healthier world.

FLY59 Rev A