

# Burkholderia cepacia complex

All USP <60> recommended strains available in EZ-Accu Shot™ and Epower™

## USP <60> recommended strains in EZ-Accu Shot™ for Growth Promotion Testing

- Quantitative microorganisms deliver 10-100 CFU per 0.1 ml
- Instant dissolve, no pre-incubation
- Traceability to reference culture ensures authenticity
- Easy & economical refrigerated storage

### Catalog Number and Strain

0488A	<i>Burkholderia cepacia</i> derived from ATCC® 25416™
01269A	<i>Burkholderia cenocepacia</i> derived from ATCC® BAA-245™
01270A	<i>Burkholderia multivorans</i> derived from ATCC® BAA-247™
0484A	<i>Pseudomonas aeruginosa</i> derived from ATCC® 9027™
0485A	<i>Staphylococcus aureus</i> subsp. <i>aureus</i> derived from ATCC® 6538™
7009	Bcc Select

## USP <60> recommended strains in Epower™ for Suitability Testing

- Quantitative microorganisms deliver 1.0-9.9E+07 CFU per pellet
- Ready-to-use format to save you time and money
- Traceability to reference culture ensures authenticity
- Easy & economical refrigerated storage

### Catalog Number and Strain

0488E7	<i>Burkholderia cepacia</i> derived from ATCC® 25416™
01269E7	<i>Burkholderia cenocepacia</i> derived from ATCC® BAA-245™
01270E7	<i>Burkholderia multivorans</i> derived from ATCC® BAA-247™
0484E7	<i>Pseudomonas aeruginosa</i> derived from ATCC® 9027™
0485E7	<i>Staphylococcus aureus</i> subsp. <i>aureus</i> derived from ATCC® 6538™

Required testing as of December 1, 2019

Visit [Microbiologics.com](http://Microbiologics.com) to place your order.

ATCC Licensed  
Derivative

\*Look for the ATCC Licensed Derivative® Emblem for all products derived from ATCC® cultures. The ATCC Licensed Derivative Emblem, the ATCC Licenses Derivative Word Mark, and the ATCC Catalog Marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

We Create Confidence in Science. Visit our website to view all strains and formats.

**microbiologics.com**

FLY133 Rev E  
2021.MAR.25

Contact us to learn more:  
320.253.7400  
info@microbiologics.com

