


# PARASITE SUSPENSIONS

CLINICAL • PROFICIENCY • EDUCATION

## ILUSTRIRANE UPUTE


1

Temeljito promiješajte suspenziju tako da ju snažno protresete ili intenzivno vorteksirate.



2


Pustite da suspenzija odleži od 5 do 10 minuta prije nego što ju primijenite u postupku testiranja.



3

Pažljivo umetnite pipetu na dno bočice i kapnite jednu kap kako biste započeli s testiranjem.

Pročitajte standardne operacijske procedure u laboratoriju, opis bojanja ili upute proizvođača za obavljanje testiranja.



 **Microbiologics®**

A safer, healthier world.