

| Demonstrates | Microbiologics Product Formats | | | | | |
|--|--------------------------------|------------|---------|------------------------|---------------------|----------------------|
| | KWIK-STIK™ | LYFO DISK™ | Epower™ | Low Inoculum Products* | Microbiology Slides | Parasite Suspensions |
| Proper Streaking Method | • | • | | | | |
| Microorganism Identification | • | • | • | | • | • |
| Media QC | • | • | | | | |
| Common Biochemical Test Methods (e.g. Catalase Oxidase MR-VP) | • | • | | | | |
| Antimicrobial Susceptibility Testing | • | • | | | | |
| Colony Morphology | • | • | • | • | | |
| Plating Technique (e.g. Membrane Filtration, Pour Plate, Spread Plate) | • | • | • | • | | |
| Growth Promotion Testing | • | • | | • | | |
| Limits Testing | • | • | | • | | |
| Proper Staining Technique | • | • | | | • | • |
| Identification of an Unknown | • | • | | | | • |

*Not intended for clinical use.



LEADING PROVIDER OF BIOLOGICAL CONTROLS FOR EDUCATION AND RESEARCH



More Than 700+ Strains in Many User-Friendly Formats



QUALITATIVE QC MICROORGANISMS

KWIKSTIK™

- User-friendly design with a lyophilized microorganism pellet, reservoir of hydrating fluid and inoculating swab – all in one device
- Over 700 strains available
- Available in packs of 2 or 6



LYFO·DISK™

- Easy-to-use: simply re-hydrate the pellet and inoculate
- Over 700 strains available
- Available in vials of 6 pellets



SLIDES AND PARASITES

QC MICROBIOLOGY SLIDES

- Contain droplets of an air-dried, fixed preserved organism, or a smear containing a designated organism population
- Wide variety of controls for instant access to uncommon specimens
- Ready-to-use out of the box of 10 slides



PARASITE SUSPENSIONS

- Perfect tool to demonstrate how to locate and identify parasite suspensions
- Over 10 parasites available including Strongyloides stercoralis and Necator americanus
- Ready-to-use vial contains 1 ml suspension
- Room temperature storage



QUANTITATIVE QC MICROORGANISMS*

e·power

- Concentrations ranging from 10² to 10⁸ CFU per pellet
- Over 76 strains available
- Ideal for schools performing food testing and water testing
- Extremely versatile
 - Can be combined for a mixed microorganism population
 - Easily manipulated to deliver a wide variety of concentrations
- Vial of 10 lyophilized pellets



LOW INOCULUM PRODUCTS*

Each product delivers 10-100 CFU per 0.1 ml

ez·accu shot™

GROWTH PROMOTION TESTING

- Instant dissolve pellets
- 5 vials of a single enumerated microorganism (1 lyophilized pellet per vial) and 5 vials of hydrating fluid (1.2 ml in each vial)
- 50 tests per kit



ez·CFU™

- 2 vials of a single enumerated microorganism (10 lyophilized pellets per vial) and 10 vials of hydrating fluid (2 ml in each vial)
- Over 900 tests per kit



ez·CFU™ One Step

- 2 vials of a single enumerated microorganism (10 lyophilized pellets per vial) and 10 vials of hydrating fluid (2 ml in each vial)
- 190 tests per kit



*Not intended for clinical use.