

EZ-ACCU SHOT™ QUALIFICATION PROGRAM

Seamless, stress-free qualification of EZ-Accu Shot for Growth Promotion Testing in your laboratory

Manual preparation of compendial strains for media challenge testing is labor intensive and presents risks in QC testing. With EZ-Accu Shot test-ready organisms from Microbiologics, you can be certain your testing and documentation will satisfy regulatory requirements and ensure QC compliance. Our one-of-a-kind EZ-Accu Shot Qualification Program is designed to simplify and guide you through the process of qualifying EZ-Accu Shot for use in the Growth Promotion Test or other media challenge testing.



How the program works:

1. Get started by completing a Qualification Program Request Form provided by your Microbiologics Sales Manager.
2. We will follow up with you and prepare a customized quote including timeline, strain requirements and pricing.
3. Your EZ-Accu Shot Qualification products are shipped to your laboratory. You will receive three individually packaged EZ-Accu Shot pellets (with hydration fluid) from three different lots of each of the five USP compendial strains and *Escherichia coli* derived from ATCC® 8739™* (18 total pellets).
4. You receive our complete guide for qualifying of EZ-Accu Shot organisms, outlining and documenting each step of the process from planning to testing and reporting.
5. Our team of Technical Support Specialists is ready to assist you each step of the way!



AVAILABLE IN THESE COMMON STRAINS AND MORE!



- **0392A** *Aspergillus brasiliensis*
derived from ATCC® 16404™*
- **0392A** *Burkholderia cepacia*
derived from ATCC® 25416™*
- **0998A** *Bacillus cereus*
derived from ATCC® 10876™*
- **0486A** *Bacillus subtilis* subsp. *spizizenii*
derived from ATCC® 6633™*
- **01092A** *Bifidobacterium animalis* subsp. *animalis*
derived from ATCC® 25527™*
- **0443A** *Candida albicans*
derived from ATCC® 10231™*
- **0487A** *Clostridium sporogenes*
derived from ATCC® 11437™*
- **0317A** *Clostridium sporogenes*
derived from ATCC® 19404™*
- **0306A** *Enterobacter aerogenes*
derived from ATCC® 13048™*
- **0366A** *Enterococcus faecalis*
derived from ATCC® 29212™*
- **0483A** *Escherichia coli*
derived from ATCC® 8739™*
- **0872A** *Geobacillus stearothermophilus*
derived from ATCC® 12980™*
- **0871A** *Geobacillus stearothermophilus*
derived from ATCC® 7953™*
- **0688A** *Kocuria rhizophila*
derived from ATCC® 9341™*
- **01090A** *Lactobacillus rhamnosus*
derived from ATCC® 53103™*
- **0419A** *Propionibacterium acnes*
derived from ATCC® 11827™*
- **0484A** *Pseudomonas aeruginosa*
derived from ATCC® 9027™*
- **0699A** *Saccharomyces cerevisiae*
derived from ATCC® 9763™*
- **0902A** *Salmonella enterica* subsp. *enterica*
serovar Choleraesuis derived from
ATCC® 10708™*
- **0421A** *Salmonella enterica* subsp. *enterica*
serovar Typhimurium derived from
ATCC® 13311™*
- **0363A** *Salmonella enterica* subsp. *enterica*
serovar Typhimurium derived from
ATCC® 14028™*
- **0360A** *Staphylococcus aureus* subsp. *aureus*
derived from ATCC® 25923™*
- **0827A** *Staphylococcus aureus* subsp. *aureus*
derived from ATCC® 6538P™*
- **0485A** *Staphylococcus aureus* subsp. *aureus*
derived from ATCC® 6538™*
- **0371A** *Staphylococcus epidermidis*
derived from ATCC® 12228™*

Why EZ-Accu Shot™ for Growth Promotion Testing?

- Traceable, test-ready organism preparations deliver 10-100 CFU per 0.1 ml as required by USP/Ph. Eur./JP
- Instant dissolve pellet reduces preparation time
- No dilutions required for simplified test procedure
- Up to 8 hours stability allows for ultimate flexibility
- 50 tests per kit, ideal quantity for mid to lower volume labs
- Refrigerated storage is easy and economical
- Online Certificate of Analysis provides detailed strain information
- Technical Support experts available for guidance



Contact us to get started

Phone: 320.253.1640 • Email: info@microbiologics.com



*Look for the ATCC Licensed Derivative® Emblem for products derived from ATCC® cultures. The ATCC Licensed Derivative Emblem, the ATCC Licensed Derivative word mark and the ATCC catalog marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

www.microbiologics.com

Contact us to learn more
320.253.1640
800.599.BUGS(2847)

 **Microbiologics®**
A safer, healthier world.